

### Table of contents

mediven®3
Active life         mJ-1       8         Select       12         Vitality       14         Travel       16         Karesse       17         Fashion Colours & Swarovski® Crystals       18         Elegance       20         Comfort       21         Active       22         For Men       23         Plus       24         Forte       25
Special needs Ulcer
Standard Harmony Glove and Gauntlet
Custom made         550/Esprit Glove, gauntlet, and sleeve
medi Accessoriesmedi Textile Care39medi Skin Care40medi Lymphpads40medi Kit Lymphedema41medi Phlebology Kit41medi Gloves41medi Fix41medi Butlers42
CircAid®. 44  Juxtacures® 46  Juxtalite® 47  Juxtalite® HD 47  CircAid® Reduction Kit 48  Juxtafit® Lymphedema Essentials 49  Juxtafit® Lymphedema Premium 50  Juxtafit® Lymphedema Sleeve 51  Juxtafit® Glove 52
CircAid® Accessories  Juxtafit® Essentials Gauntlet
medi e-chon







# We have the direct connection

To communicate. To co-operate. To be there for each other. The medi service team respects and values this requirement seriously.

We are always here for you. This is your medi service team.

#### **Customer Service**

for consultancy, ordering, problem solutions

Tel: 1800 361-3153 Fax: 1888 583-6827

You can visit our website:

contact us and place orders by email:

service@medicanada.ca

Place orders online 24h/7 on our online boutique:

eshop.medicanada.ca

www.medicanada.ca



### mediven® features



#### Clima Comfort\*

- Wicks moisture away from the skin's surface quickly
- Guarantees a high level of breathability
- \* Now scientifically tested (Institut Francais du textile et de l'habillement (IFTH), Clima Comfort Test, Villeneuve-d'Ascq, 2010.)



#### **Integrated Clima Fresh function**

- Antibacterial effect
- Thus inhibits odour-forming bacteria



#### Perfect Fit

- A very precise fit
- Medically effective pressure gradient
- Optimal care conditions in standard and made-to-measure version



A wide range of variants and special models



### Anatomical oblique borders as standard (distal and proximal)

- Unimpeded lymphatic drainage without circular constriction
- · Optimum fit and better hold

#### Colours for more quality of life

- Wide range of standard colours readily available
- Trend colours available within 15 working days

#### Easy to care for

 Suitable for the delicates program of washer and drying machines







All mediven stockings are latex-free.

# mediven® Standard colours overview



STANDARD COLOURS	MEDIVEN KARESSE	MEDIVEN ELEGANCE	MEDIVEN COMFORT	MEDIVEN PLUS	MEDIVEN FORTE	MEDIVEN FOR MEN	MEDIVEN ACTIVE
Cashmere Poudre <sup>1</sup>	1						
Caramel Ficelle <sup>2</sup>	2						
Noir		•	•	•	•	•	
Anthracite Mouette³	3		*	*	*		
Navy			*	*	*		
Brown							
White		0		*			0
COLOURS AVAILABLE AS SPECIAL ORDERS	MEDIVEN KARESSE	MEDIVEN ELEGANCE	MEDIVEN COMFORT	MEDIVEN PLUS	MEDIVEN FORTE	MEDIVEN FOR MEN	MEDIVEN ACTIVE
Rosé		*	*	*	*		
Bronze		*	*	*	*		
Beige		*	*	*	*		
Magenta		*		*			

<sup>\*</sup>Special order: delivery period up to 15 working days. Non-returnable.

### mediven®

#### SUPPORT STOCKINGS FOR PREVENTION

### **m** J-1

Class: 15-20 mmHg Knits: Metropole, City and Cato

Colours: Variety of colours depending on knit Styles: Calf, Thigh and Pantyhose

For Women and Men

**Our Ultra-Sheer** 

For Women and Men Travel

Class: 18-20 mmHg Style: Calf length Colours: Black and Anthracite for men / Women: Black and Caramel

For Women and Men Vitality, Select

Foot: Closed Class: 15-20 mmHg, CL I

Style: Standard in calf

#### FOR CASES OF HEAVY LEGS. FATIGUE WITHOUT SWELLING **Karesse**

Class: CLI Foot: Closed

Styles: Calf, Thigh and Pantyhose

Colours: Black, Powder, Ficelle and Mouette

#### FOR CONTROL OF MILD EDEMA

#1 for Women **Elegance** 

Classes: CLI - CLII, Open or closed foot

Sizes: Standard and Custom

Very wide range of styles, patterns, colours and Swarovski® crystals

#### Comfort Rated number 1 for comfort

Classes: CL I, Closed Toe Colours: Cashmere, Caramel and Black

Styles: Calf, Thigh and Pantyhose Trend colours available

#### FOR CONTROL OF MODERATE EDEMA

For the active person **Active** 

Classes: CL I - CL II, Standard sizes only

Style: Ribbed calf sock with a sporty look Colours: Navy, Black, Brown, Anthracite & White

For men For the everyday look

Classes: CLI-CLII Foot: Anatomical Fit

Colours: Navy, Black, Brown and Anthracite Style: Calf only, Modern flat ribbed look

Plus When more control is required

Classes: CL I - CL II - CL III Foot: Open or closed

Styles: Unisex Wide range of styles, Our all rounder

#### FOR SEVERE EDEMA CONTROL STARTING EARLY IN THE DAY

**Forte** 

Classs: CL II - CL III

Custom Made Flat Knits: Mondi, Cosy, 550

Classes: CLI - CLII - CLIII - CLIV

























The ultimate control in our circular knits

Foot: Open or closed

Style: Robust knit with antibacterial treatment throughout the stocking

Colour: Variety of colours and styles

The Ultimate choice!

Foot: Open or closed

Styles: Variety of flat knits in a vast choice of styles and options





#### Colours







Pearl

### **MJ-1** Metropole

### Class: 15-20 mmHg

#### Closed toe

Styles: Calf, Thigh with silicone band and Panty

Size: S, M et L

Lengths: L1 et L2

The Ultra sheer fully measured mJ-1 is the knowledgeable and discreet choice due to it's aesthetics and efficiency.

woman: MJ-1 City in Midnight | man: Mj-1 City in Jet Black







Colours







Gradual introduction of new colours: Oxid (replaces Basalt) Nerea (replaces Midnight) Colours





Gradual introduction of new colours: Nerea (replaces Midnight)

### **MJ-1** City

#### Class: 15-20 mmHg

#### Closed toe

Style: Calf

Sizes: S, M, L

Lengths: L1 and L2

The City support sock was made for Men, but women love it too, for its opaque appearance. The microfiber knit is soft and easy to wear.

### **MJ-1** Cato

### Class: 15-20 mmHg

#### Closed toe

Style: Calf

Sizes: S, M, L, L+

Lengths: L1 and L2

The Cato support sock for men and women, has a textured finish and is both opaque and very robuste.









### Why they're different.









#### Luxurious fiber blend

Knit integrally into every Rejuva garment.

#### **Fashion patterns**

Perfect for work or dressy occasions.



### Reciprocated heel & toe

Reduces friction for added comfort.

## New!



## mediven® select

Premium compression for every walk of life!

#### **Features:**

- Casual pattern for work and leisure.
- Plush fabric, constructed of high-quality spandex and added Lanolin, for unparalleled softness.
- Provides accurate, consistent, long-lasting medically graduated compression.
- Soft, wide silicone top band for a comfortable, stay-put fit.
- Distinct left and right foot designs for improved fit.



### mediven select

#### Casual, patterned sock for work or play

#### Indications:

- · Mild to moderate spider veins.
- · Slight to moderate varicose veins.
- · Mild edema.
- · Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

### 

#### material composition

77% polyamide, 23% elastane

#### washing instructions



Machine wash and dry.



select knee-high sock in black



select pattern showing the wide top-band



classic pattern showing anatomically shaped foot with right and left markings

Circumference	II	III	IV	V	VI	VII		
CALF	30 – 39 cm 11¾" – 15¼"	32.5 – 42 cm 12¾" – 16½"	35 – 45 cm 13¾" – 17¾"	37 – 48.5 cm 14½"– 19"	39.5 – 51.5 cm 15½" – 20¼"	41.5 – 54.5 cm 16¼" – 21½"		
EXTRA-WIDE CALF	N/A	42 - 49.5 cm 16½" - 19½"	45 - 53 cm 17¾" x 20¾"	48 - 56 cm 19" x 22"	51.5 - 59 cm 20¼" - 23¼"	54.5 - 62 cm 21½" - 24½"		
ANKLE	20.5 – 22.5 cm 8"– 8¾"	23 – 25 cm 9"– 9¾"	25.5 – 27.5 cm 10"– 10¾"	28 – 30 cm 11"– 11¾"	30.5 – 33 cm 12"– 13"	33.5 – 36 cm 13¼"– 14¼"		
LEG LENGTH	calf > standard: <47.5 cm (<18") tall: <47.5 cm (>18")							

<sup>\*</sup> Extra-wide calf available in black only

## New!



# mediven® vitality

Premium compression for every walk of life!

#### Features:

- Luxurious patterned knee-high socks for men and women.
- Herringbone-ribbed design delivers a classic look.
- Enhanced softness for all-day comfort.
- · Medi-grip top-band provides seamless security without irritation.
- Style and feel maintains top-tier medical efficacy you expect from medi.

#### **Benefits:**

- Clima Fresh technology inhibits bacteria growth for clean, refreshed legs.
- Clima Comfort blended with Lanolin keeps legs comfortably dry all day.
- Available in standard and petite lengths.



### mediven vitality



#### Your on-the-go compression stocking

#### **Indications:**

- · Minor to moderate spider veins.
- · Slight to moderate varicose veins.
- · Mild to moderate edema.
- Tired, achy, heavy legs.
- Helps prevent deep vein thrombosis (DVT).

#### colours

CHARCOAL

EBONY

compression levels

I II III IV V VI VII

15-20, 20-30 mmHg

#### lengths

width

petite and standard

standard

#### material composition

77% polyamide, 23% elastane

#### washing instructions












Machine wash and dry.



Vitality knee-high sock in charcoal



Herringbone-ribbed pattern detail



mediGrip top-band technology

Circumference	I	II	III	IV	V	VI	VII	
2 CALF	28 – 35.5 cm	30 – 39 cm	32.5 - 42  cm	35 – 45 cm	37 – 48.5 cm	39.5 – 51.5 cm	41.5 – 54.5 cm	
	11"– 14"	11³/4"– 15¹/4"	$12^{3}/_{4}" - 16^{1}/_{2}"$	13³/₄"– 17³/₄"	14¹/2"– 19"	15 <sup>1</sup> /2" – 20 <sup>1</sup> /4"	16 <sup>1</sup> / <sub>4</sub> " – 21 <sup>1</sup> / <sub>2</sub> "	
1 ANKLE	18 – 20 cm	20.5 – 22.5 cm	23 – 25 cm	25.5 – 27.5 cm	28 – 30 cm	30.5 – 33 cm	33.5 – 36 cm	
	7" – 7³/4"	8"– 8 <sup>3</sup> / <sub>4</sub> "	9"– 9³/4"	10"– 10³/4"	11"– 11³/4"	12"– 13"	13 <sup>1</sup> / <sub>4</sub> "– 14 <sup>1</sup> / <sub>4</sub> "	
LEG LENGTH								





For him: Ribbed finish









### mediven® Travel For Him and Her

### Class: 18-20 mmHg

Closed toe

Style: Calf

Sizes: S, M, L, XL, XXL

Lengths: Petite and Regular

Who doesn't travel nowadays? People spend a lot of time travelling in small confined spaces, whether in a bus, train, plane or car. mediven® Travel was developed by a team of physicians, technicians and end users which guarantees a highly effective and user-friendly product. mediven® Travel is the only clinically tested prophylaxis travel garment on the market today!





NEW COLOURS

NOW AVAILABLE!

#### Colours









Powder available soon

### mediven® Karesse

#### Class: I

#### Closed toe

Styles: Knee, thigh with silicone lace and panty

Sizes: I - V

Lengths: Petite and Regular

Our trendy ultra-sheer stocking with its push up panty will slim and lift your figure. The concept behind mediven® Seduction is it's innovative knit: a unique, soft, breathable knit with no unattractive seams, a nude toe and a discreetly reinforced heel. The knit is extremely soft.



Say
hello
to our
new
trend
colours



Compression everyone wants to wear.

- New fashionable trend colours
- New motifs and crystals
- Wide variety of functional options
- Effective, high-quality products
- Suitable solutions for virtually all phlebological and lymphological indications

# mediven® 2020 trend colours overview

TREND COLOURS	MEDIVEN ELEGANCE	MEDIVEN COMFORT	MEDIVEN PLUS	MEDIVEN FORTE
Cherry-red	*	*	*	*
Blue-jeans	*	*	*	*
Violet	*	*	*	*
Grey	*	*	*	*

### Crystals from Swarovski®\*\*


### option available for all mediven main qualities



Special order: Swarovski® Crystals\*\* can be added to all mediven® flat knit products and to the following circular knit products: mediven elegance, mediven comfort, mediven plus and mediven forte.

<sup>\*</sup>Special order: delivery period up to 15 working days. Non-returnable.

<sup>\*\*</sup>Swarovski® is a registered trademark of Swarovski AG.



mediven® Elegance

Classes: | - ||

Closed toe (Open Toe available in Cashmere Only)

**Styles:** Calf, Calf SB, Thigh, Extra-wide Thigh/SB, Panty and Maxi queen size, Maternity.

Sizes: I - VII

Lengths: Petite and Regular

The new innovative hybrid knit technology combines transparency and softness into one product: mediven Elegance. Available in the widest range of styles. Colours, motifs and crystals

**Maximum sheerness** 

» Makes the product easier to wear


Hybrid knit technology combines softness and transparency

» Makes the stocking easier to apply and fit like a second skin

STANDARD

CUSTOM MADE





### mediven® Comfort

#### Class: |

#### Closed toe only

Styles: Calf, Thigh with silicone band and Panty

Sizes: I - VII


Lengths: Petite and Regular

Rated Number 1 by the Hohenstein Institute for it's comfort and it's lanolin based moisturising properties. This stocking makes you forget that you are wearing a medical grade compression stocking. Only the best for you! mediven® comfort has a unique ultra soft feel while providing a perfect fit. Double-knitted comfort sole for a more pleasant walking sensation.

STANDARD

CUSTOM MADE





Colours











### mediven® Active

#### Classes: I - II

Closed Toe

Style: Calf

Sizes: I - VII

Lengths: Petite and Regular

Turn men into faithful customers with new mediven® Active. Your female customers will love it too. mediven® Active is a robust compression garment with a ribbed design that withstands unrelenting wear and tear - perfect for the physically active. It is also impressive for its medical effect, comfort and softness. Soft heel, Re-inforced sole and robust antibacterial treatment.





Colours









### mediven® For Men



Anatomical foot for better comfort

#### Classes: | - ||

Closed toe

Style: Calf

Sizes: I - VII

Lengths: Petite and Regular

Made for him but loved by her. A modern flat ribbed look and a unique anatomical foot, offers great comfort while providing hydrating properties. Our patented technologies Clima-Comfort (moisture wicking knit) and Clima Fresh (antibacterial process) have been optimised in our mediven® for Men.

























White

STANDARD

CUSTOM MADE Trend colours available

### mediven® Plus

#### Classes: | - || - |||

#### **Closed or Open toe**

Sizes: I - VII

Lengths: Petite and Regular

Wide range of styles, 3 standard colours (beige and black), in ready-wear and custom-made. Other colours are in special order.

Mediven® Plus is the best comfort alternative when a firmer knit is required. It provides patients with fashionable graduated compression therapy in a wide range of ready-made styles. The ideal solution and problem solver for fitting issues.

**Special Options:** Calf with silicone band available as a standard product.

Extra-wide calf available for classes 1 and 2.





Standard colours

Colours available as special order

Caramel Cashmere

Black

Anthracite

Navy

Rosé

Bronze

Beige

Trend colours available

### mediven® Forte

Classes: || - |||

STANDARD + CUSTOM MADE

**Closed or Open toe** 

**Styles:** Calf, Thigh with or without silicone band Thigh with waist attachment

Sizes: I - VII and custom Made

Lengths: Petite and Regular

Wide variety of styles. Available as custom-made with special options. Mediven® Forte is our strongest circular-knit compression stocking, guaranteed wrinkle free.

It is used for the treatment of pronounced venous disease, in the early stages of lymphedema, healing ulcers, severe disorders of the venous function and after thrombosis.













### mediven® Ulcer

#### Class: 2 x 20 mmHg

Closed and open Toe Style: Calf


Sizes: I - VII

Lengths: Petite and Regular

The special feature of mediven® Ulcer stocking is the silver that is incorporated into the yarn. Silver, even in low concentration, is an effective barrier against a wide range of bacteria and also promotes healing, making this stocking an excellent choice when treating ulcers.

The mediven® Ulcer Kit consists of 2 compressive components. One pair of mediven® Ulcer stockings with an ankle pressure of 20 mmHg designed for continuous around-the-clock (20 mmHg) compression and mediven® Ulcer Plus for daytime wear (with ankle pressure of 20 mmHg for combined pressure of 40 mmHg). The use of both elements guarantees maintenance of prescribed pressure at all times, whether at rest or active.

Circumferen	ce	1	II	III	IV	V	VI	VII
Ankle	B	18 - 20 cm	20 - 22 cm	22 - 24 cm	24 - 26 cm	26 - 28 cm	28 - 30 cm	30 - 32 cm
Calf	0	28 - 34 cm	30 - 37 cm	33 - 40 cm	35 - 43 cm	37 - 46 cm	39 - 49 cm	41 - 51 cm
Lenght		Petite	Regular					
Calf	0	34 - 38 cm	39 -44 cm					



	Le	Colour	tagging					
Size	Ankle	Calf	Thigh	Length	Ankle	Thigh		
	D	G	G	0				
Knee								
XS	17 - 19	26 - 32	n/a	35 - 46	orange	n/a		
S	20 - 22	30 - 38	n/a	35 - 46	green	n/a		
M	23 - 25	34 - 44	n/a	35 - 46	violet	n/a		
L	26 - 28	38 - 48	n/a	35 - 46	blue	n/a		
L-X	26 - 28	44 - 55	n/a	35 - 46	grey	n/a		
XL	29 - 31	44 - 55	n/a	35 - 46	red	n/a		
			Thigh					
XS	17 - 19	26 - 36	40 - 50	60 - 84	orange	orange		
S	20 - 22	30 - 42	48 - 60	65 - 89	green	green		
S-X	20 - 22	30 - 42	56 - 70	65 - 89	green	grey		
M	23 - 25	34 - 46	56 - 70	65 - 89	violet	violet		
М-Х	23 - 25	34 - 46	64 - 80	65 - 89	violet	grey		
L	26 - 28	38 - 51	64 - 80	65 - 89	blue	blue		
L-X	26 - 28	38 - 51	72 - 90	65 - 89	blue	grey		
XL	29 - 31	42 - 55	72 - 90	65 - 89	red	red		
XL-X	29 - 31	42 - 55	80 - 100	65 - 89	red	grey		
		Thigh	n with hip f	ixation				
S	20 - 22	30 - 42	48 - 60	65 - 89	green	green		
М	23 - 25	34 - 46	56 - 70	65 - 89	violet	violet		
L	26 - 28	38 - 51	64 - 80	65 - 89	blue	blue		
L-X	26 - 28	38 - 51	72 - 90	65 - 89	blue	grey		











### mediven® Thrombexin 18 Anti embolism

### Class: 18 mmHg


#### **Closed and Open Sole**

Styles: Calf, Thigh with silicone band, Thigh with waist attachment

With a pressure of 18 mmHg mediven® thrombexin is effective alternative to traditional postoperative compression bandages or for short-term treatment.

Simple, fast, safe and effective. The system of colour-coded sizes can quickly and easily validate the size of the garment. Ease of inspection and treatment, thanks to an inspection opening. This medical stocking is suitable for thigh circumferences of up to 100 cm.

Widest sizing on the market, best fitting anti embolism stocking.







# Lymphedema

### Circular knit

#### Harmony For light to medium control of lymphedema

**Knit:** Circular textured knit for excellent adherence to the skin, provides extra comfort at elbow

Classes: CLI-CLII

Styles: Gauntlet, Glove, Sleeve and Combination sleeve with or without silicone,

Sleeve and Combination sleeve with shoulder attachment

Sizes: Standard, offered in regular and X-wide

Colours: The standard colour is Caramel. Certain styles are also available in Sand

### **Standard**



### **Custom Flat Knit with Seam**

#### **Esprit** For medium to strong control of lymphedema

Knit: Flat knit with seam with possibility of colours, 3 patterns

and Swarovski Crystals

Classes: CLI-CLII-CLIII

Mondi

Sizes: Custom Made, for upper extremity

and Swarovski Crystals

Classes: CLI-CLII-CLIII

Colours: See custom order form for full colours

Knit: Flat knit with seam with possibility of colours

### **Custom Made**





### Cosy For medium to severe control of lymphedema

For medium to severe control of lymphedema

Knit: Flat knit with seam with possibility of colours, 3 patterns

and Swarovski Crystals

Classes: CLI - CLII - CLIII

Sizes: Custom Made only, for lower extremity

**Sizes:** Custom Made only, for lower extremity **Colours:** See custom order form for full colours

Colours: See custom order form for full colours

#### 550 Stockings/Sleeve For optimal Control of severe lymphedema

Knit: Flat knit with seam with possibility of colours, 3 patterns

and Swarovski Crystals


Classes: CL I - CL II - CL IV (leg only)

Sizes: Custom made for upper and lower extremity

Colours: See custom order form for full colours









### mediven® Harmony

#### Classes: | - ||

Knit: Circular textured knit for excellent adherence to the skin, provides extra comfort at elbow

Styles: Gauntlet, Glove, Sleeve, Combined sleeve and Sleeve with shoulder attachment

Sizes: Standard, offered in Regular and Extra-Wide

\* Colours: offered in 2 colours in the following styles: regular sleeve with silicone band and combined sleeve with silicone band. Availability of sizes may vary

mediven® Harmony is the latest range of medi sleeves. It raises the standards in circular knitting with a perfect balance between efficiency and comfort. It is knitted with high quality standards to ensure accuracy and graduated compression. mediven® Harmony is a unique fabric that is soft to the touch, clinically effective and easy to apply. It features a new design that incorporates a double ellipse - one in the crease of the elbow and the other on the tip of the elbow - for maximum comfort to sensitive areas.

Harmony is used in lymphedema stage 0-1 controled.



### mediven harmony glove & gauntlet

New seamless version **now** available!

Matched to the mediven arm garment, medi offers the appropriate hand pieces and gloves. The mediven harmony compression glove and gauntlet offer patients maximum wearing comfort and make the therapy simple and pleasant.



### **Optimal wearing comfort**

Due to its soft and breathable material, the glove and gauntlet are gentle on the skin and cause no skin irritation.

As with the arm garment, Lanolin has also been added to the glove and gauntlet to provide the benefits of natural moisturizing for a soft and smooth feel.



### Therapeutic transition zone

The sleek and soft top border ensures a discreet transition to the arm garment, complementing it perfectly.







#### Precise compression for sensitive points

NEW

The seamless construction of the new harmony glove and gauntlet reduces pressure points that seams can cause while allowing more freedom of movement.

### mediven harmony® armsleeve The proven benefits

### **Freedom of Movement**

The innovative fabric within the double ellipses eases the feeling of pressure on the point and the crease of the elbow, making flexion and extension of the elbow comfortable and unhindered.

The ellipse in the elbow crease encourages optimum lymphatic drainage. Fabric wrinkling is reduced to a minimum to prevent rubbing and skin soreness.

# Proven therapeutic compression

The innovative fabric is guaranteed to be effective. Especially adapted for post-operative and post-traumatic treatment and even suitable for softer connective tissue; mediven harmony combines high wall stability with a gentle touch from the supple knit.

### Luxuriously soft

The surprisingly soft feeling is the result of new product developments. In addition to the double ellipse the innovative fabric structure stimulates the tissue and exerts a pleasant massage effect.

Lanolin has been added to the armsleeve to provide the benefits of natural moisturising for a soft and smooth feel.

### Security and comfort

The narrow beaded silicone topband is discreet and effective. The topband is knitted to the armsleeve with a flat seam. The silicone pattern offers security with breathability. Armsleeves are also available with a knitted border.



mediven harmony: ellipse in the elbow crease



mediven harmony: ellipse on the point of the elbow



mediven harmony topband



Black



Sand












Trend colours available

CUSTOM MADE





### mediven® Esprit / 550 honeycomb knit

**Custom Made Esprit / 550** 

Classes: | - || - |||

Knit: Custom made flat knit with seam, possibility of colours and three designs

Styles: CUSTOM MADE: Unlimited possibilities

Sizes: Custom Made

Mediven® Esprit is used for the treatment of lymphedema stages 1 and 2, in the cases of post-operative and post-traumatic edema, and after burns. Knit provides micro massage.

#### **Custom Made**














And other options ...



CUSTOM MADE

### mediven® Mondi honeycomb knit

Classes: | - || - |||

Knit: Custom flat knit

Styles: Calf, Mid-thigh, Thigh, Panty, One leg panty, Cyclist shorts, Capri and Leggings

mediven® Mondi provides optimal fit and therapeutic action. Possibility of incorporating pads, Silicone band, Velcro closures and zippers.

Optional Individual knitting marks facilitate donning





### mediven® cosy

Flexible. Supple. Freedom of movement.

#### Suitable for:

- Soft to firm connective tissue
- Moderate deepened skin wrinkles / calibre jumps
- Associated weight disorders
- · Lipoedema patients
- Travelling / sports

#### Indications

- · Lipoedema stages 1, 2, 3
- Primary and secondary lymphoedema stages 1, 2
- Severe venous indications with oedema tendency
- Scar treatment
- Post-operative and post-traumatic oedema
- · Mixed forms such as phlebo-lymphoedema, lipolymphoedema

Design-Elements\*



Flower



Lilac



Compression

CCL 1, 2, 3

Oedema stage


**Product benefits** 

Clima Fresh



#### **Care instructions**

Easy to care for <u>@</u> ⊗ 🔉 🗷 🖸 Suitable for the delicates programs of washing machines and dryers.





Trend colours available



<sup>\*</sup>The Design-Elements can be freely combined with the standard and trend colours (no extra charge).

\*\*Swarovski® is a registered trademark of Swarovski AG. Swarovski® crystals can be freely combined with the standard and trend colours.

\*\*\*Available as seamless toe cap or individual toe cap.





















Trend colours available

CUSTOM MADE

### mediven® 550 honeycomb knit

### **Classes:** | - || - ||| - |V

Knit: Custom flat knit with seam, possibility of colours and three patterns

Wide range of options: Custom made for upper and lower extremity

Of All our flat knit garments, mediven® 550 is the most classic and versatile. Ideal for lymphedema, it offers a micro massage. It presents a wide range of options and is available in several colours and patterns. Possibility of incorporating pads, Silicone band, Velcro closures and zippers.

550 is used in lymphedema stage 1 - and over.

### Design-Elements







#### Fashion-Elements\*




















# mediven® mondi\*\*, cosy & 550 toe cap

#### **Product features**

- Toe cap with Clima Fresh: reliably prevents odour and colonisation by bacteria
- Top and bottom principle:
  - Seams moved to the upper and lower surfaces of the toes in the direction of lymphatic drainage
  - Lateral seams give every toe stability for an effective action Each single toe made to measure anatomically for an optimum fit
- Optional integration of medi lymphpads for additional pressure in lymphoedema of the dorsum of the foot
- \*\*Fashion-elements not available for mediven mondi



\*The Design-Elements can be freely combined with the standard and trend colours (no extra charge).





M = Made-to-measure

Trend colours available



Seamless toe cap

#### Seamless toe cap

- More reliable measurements (no longer necessary to calculate the overlap)
- Care to suit the indication for optimum lymphatic drainage (no circular seams)
- Very comfortable to wear, without circular seams
- Outstanding non-slip grip
- Available in variants AB, AB1, AD, AF, AG, AT, ATH, ATU, ATE/ATEH
- Different compression classes in toe cap and leg garment available



Individual toe cap

#### Individual toe cap

- Increased wall stability due to double compression garments for pronounced forms of oedema over the dorsum of the foot
- Increased treatment compliance due to quick and easy donning
- Individual care options, e.g. nightwear or combination with other product variants
- Displacement of the seam on the bottom of the foot (lateral seam as special request available)



# How to Care for your compression garments ...

### For the Hygiene of your garments

Mediven® offers a complete line of cleaning products and skin care to meet your needs.

- Wash your compression garments daily to benefit from their full effectiveness. Spacing out the washings will not extend their life.
- As an alternative to mediven® soap use a liquid detergent such as Tide, Sunlight, Arctic Power, etc ... without bleach or fabric softener.
- All other types of soap such as shampoo, dish soap, liquid hand soap and bar soap should always be avoided as they destroy the elastic fiber of your garment. Only garments without silicone can be machine dried at low temperature for up to 5 minutes.
- Never choose a laundry detergent that contains bleach or fabric softener.
- Do not use soap for delicate laundry as Zero, Woolite, Forever New, which already contain softeners.

- Turn your stocking inside out when washing, this
  will help release any residual dead skin and clean
  the silicone band. Turn stocking back to the correct
  side and hang to air dry in open air. No machine
  drying for garments with silicone band
- All mediven® seamless stockings or garments can be washed with warm water either by hand or machine washed on the gentle cycle in a washing bag.
- Washing your Flat knit garments with seams may be hand washed or machine on delicat cycle. Once washed, stretch the garment out so that it can regains its original dimensions then hang to dry.

FOR FURTHER INFORMATION, PLEASE CONTACT OUR CUSTOMER SERVICE.

### medi Textile Care

### Gentle care for hand washing your garments



medi Clean

**Tube: 150 ml** 

For the easy care of your compression garments -anytime, anywhere

The special washing liquid is specifically compatible for use with compression garments, cleaning with the utmost of care.



medi Spot ex

The clever solution for stain-free compression garments

Use before the wash, the spot removal stick medi spot ex removes shoe marks from the garment's toe or heel..

Roll: 15 ml



### medi Laundry net Gentle washing

Gentle care helps compression stockings maintain their medical efficacy.

The medi laundry net is the ideal aid for caring for compression stockings. The net structure protects the fabric against damage during the daily wash, but does not hinder the removal of dirt particles or soap residues.

### medi Skin care

Venous and lymphatic diseases are often accompanied by an imbalance in the skin. The congestion of blood and fluid in the leg can also adversely affect the skin. To enhance the success of compression therapy, three products were developed to meet specific requirements.







### medi Day Gel

### Tube: 50 ml

### medi Night Care

#### Tube: 50 ml

### medi Fresh spray

### bottle: 100 ml

Medi day, cools and freshens throughout the day.

Reduces the tension and sensitivity of the skin

Support skin regeneration

Apply on the skin a few minutes before putting on the socks

Heals and relaxes the superficial layers of the epidermal

Helps regulate water balance

Accelerates regeneration and improves skin appearance

Apply to the skin after removing compression stockings and gently massage into the skin.

Refreshes and cools

Ideal complement to medi Day gel.

Ideal for work or travel

Spray directly onto the compression stocking or the skin several times a day as required. The spray will not damage the knit of the stocking.

### medi Lymphpads

medi lymph pads for the arms stimulate the microcirculation and increase lymph outflow. Hardened tissue or areas of fibrosis are loosened to provide general support in treatment of oedema.







### mediven® Lymphpads

Cushioning sheet Lymphedema - Ribbed (2 per pack)

Cushioning sheet Lymphedema - Fine Dots (2 per pack)

Cushioning sheet Lymphedema - Big Dots (2 per pack

Buckwheat cushion shaped kidney

Flaxseed cushion shaped kidney

**Buckwheat cushion shaped Pear** Flaxseed cushion shaped Pear

Buckwheat cushion Round shape

Flaxseed cushion Round shape

Cushions



### medi Vascular Kit

This Kit contains the following items:

- 1 Pair of Medi Donning Gloves
- 1 Sample of Medi Comfort Night Cream
- 1 Sample of Medi Comfort Day Gel
- 1 Sample of Medi Clean Soap
- 1 Patient Guide



medi Lymphology Kit

This Kit contains the following items:

- 1 Pair of Medi Donning Gloves
- 1 Laundry Bag
- 1 Sample of Medi Comfort Night Cream
- 1 Sample of Medi Comfort Day Gel
- 1 Sample of Medi Clean Soap
- 1 Measuring Tape
- 1 Lymphedema Patient Guide



medi Donning gloves Available per/Pair or box of 12

Sizes: S, M and L

Rubber gloves for perfect adherence

Donning stocking made easy! Rubber gloves for perfect adherence. These gloves with a rough inner surface allows good grip.



medi Textile Gloves
Available per/Pair or box of 12

Sizes: S, M, L and XL

Fingertips coated with silicone up to the distal phalanx for a perfect adherence

These gloves with a coarse inner surface will help position the stocking on the leg. Donning stockings made easy!



medi Fix medi fix with 50 ml (unit with 6 pieces)

Adhesive lotion for securing medical compression stockings and elastic compression garments.

### medi Butler

#### Benefits for customers

Quick and simple donning of compression stockings Less effort required Wide choice for the most individual demands





### medi Export Butler

#### With long handles

The ideal donning aid for all compression stockings wearers who can comfortably don their stockings without assistance. The long handle ensures maximum stability when putting on compression stockings and pantyhose.

#### medi Butler variable handles

#### With variable handles

The donning aid for stocking wearers with considerably restricted freedom of movement.

The handle length can be adjusted between 37 and 58 cm.





### medi Big (XXL) Butler

#### For large circumferences

#### Sizes IV to VII or customized

The stirrup of the medi Big Butler has a larger diameter. This makes allowance for larger circumferences which are often associated with less freedom of movement.

### medi Arm Butler

#### With short handles

#### For sleeves

Perfectly suited for donning armsleeves or for helping a carer put compression stockings on to the person being cared for.



### medi Big Arm Butler

#### Arm-Big Butler

#### Available in 2 sizes


The medi Arm Butler eases donning of compression armsleeves. This unique arm butler hooks on to flat surfaces for better stability. Available in regular and large size.



### medi Hospital Butler

#### Help with threading

Adjustable diameter with hard paint finish coating.



Simply put on

compression stocking

### medi Butler 2 in 1

#### Assist with threading and removal

Donning and doffing aid expertly designed to get you in and out of your stockings in no time. Discreet, it fits easily in a pocket or purse and provides full autonomy anywhere you go. Suitable for open or closed toe stockings as well as arm sleeves.





### medi Butler Off

### The perfect assistant to remove her stockings!

#### Out of the stocking in seconds

The medi Butler Off helps you take off your compression stocking without tiresome bending and without having to pull hard. To take the stocking off, simply slide the medi Butler Off into the band and push down.





## CircAid® by medi

### Venous

**Juxtacures**<sup>®</sup>

Venous pathology (ulcer)/calf

Juxtacures® Lower legging

Comfort™ Leg liner (without foot)

Comfort™ Compression anklets (pair)

Available in 3 lengths with a justable circumference that can follow the evolution of the patient



Juxtalite® Standard & HD

Venous Insufficiency ulcers / Calf

Juxtalite® Lower legging C

Circaid Compressive Undersock

Juxtalite® HD Lower legging



### Lymphedema

CircAid® Reduction Kit

**Reduction Kit Lower leg** 

Reduction Kit Knee

**Reduction Kit Armsleeve** 



Juxtafit® Essentials

Juxtafit® Essential Lower leg

Comfort™ Leg liner (without foot)

Comfort™ Compression anklets (pair) with Pac Band



Juxtafit® Premium

Calf - Knee - Thigh

Lower legging available in standard. Calf, knee, thigh or combined available in custom made

Juxta t® Premium Legging

Comfort™ Silver leg liner (without foot)

Comfort™ Compression anklets (pair) with Pac Band



Juxtafit® Essentials

Sleeve

Available in custom made only

Juxtafit® Custom Armsleeve

Comfort™ Coverup



**Juxtafit®** Essentials

Glove

Available in custom made

Juxtafit® Glove

**Option:** Custom Made Dorsum Strap





### **Juxtacures®**

### Classes: 20-30, 30-40 and 40-50 mmHg

#### Warranty: 6 months / manufacturing defect

Juxtacures® is a instantly adjustable, inelastic multi-layer compression system that may be used by hospitals and clinicians as a secondary dressing for the treatment of open venous stasis ulcers. This patented bandaging alternative is easy to use, comfortable to wear and can help improve healing rates. Juxtacures® Built-In pressure System allows patients with venous stasis ulcers or mixed arterial-venous disease to set compression to an acceptable prescribed level. Verifying the correct level of applied pressure is simple with the BPS guide card, just measure, verify and adjust – you can't do that with bandages.

This product can be worn 24h/7 days a week. Juxtacures® can be modified to follow the reduction of edema.

 $^{\star}$  The difference of the value 1 and the value 2 must not exceed 30 cm.

Circumference		Measuring points	Short	Standard	Long
Length	0	Heel with knee bend	- 44 cm	44 - 49 cm	+ de 49 cm
Circumference	G	Calf	minim	num 25 cm maximum	64 cm
Circumference	erence B Ankle minimum 19 cm		num 19 cm maximum	42 cm	
Circumference	A)	Arch of the foot	minim	num n/a maximum	34 cm





- 1 CircAid black compressive undersock
- 1 Built-in Pressure guide card (BPS)

- 2 CircAid Black Compressive Undersock
- 1 Built-in Pressure guide card (BPS)



### Juxtalite® Venous Standard

### Juxtalite® Venous HD

### Classes: 20-30, 30-40 and 40-50 mmHg

#### Warranty: 6 months / manufacturing defect

Juxtalite® is best know as an alternative to elastic compression stockings for the treatment of chronic venous insufficiency, it is also an excellent product for the treatment of mild to moderate Lymphedema of the lower leg. Because of its thinner and lighter material Juxtalite® is not recommended for severe cases of lymphedema. Juxtalite® is an inelastic adjustable compression garment made of "Breath-O-Prene" and can be worn day and night. Available in 8 standard sizes and 2 leg lengths. Juxtalite® can fit most patients. This alternative permits patients to apply the prescribed compression. Juxtalite® has a warranty of 6 months.

This product can be worn 24h/7 days a week.

#### Indications: moderate edema, CVI, active & healed venous ulcers, phlebolymphedema

The new juxtalite hd (heavy duty!) is made of a thicker, more robust material, making it better suited for patients with more extreme edema.

### Size chart for juxtalite & juxtalite hd

Circumfere	ence	S	M	<b>M</b> wide calf	L	<b>L</b> wide calf	XL	XL wide calf	XXL
Ankle	В	19 - 22	22 - 27	22 - 27	25 - 32	25 - 32	30 - 37	30 - 37	35 - 42
Calf	G	26 - 36	33 - 44	44 - 54	41 - 51	51 - 61	37 - 48	48 - 61	53 - 64
Lenght									
Low knee	.ow knee								



1 Reduction Kit Calf

- 1 Reduction Kit Calf 1 Reduction Kit Knee
- 1 Reduction Kit Tigh

1 Reduction Kit Arm 1 gantelet 1 comfort sleeve

### CircAid® Reduction kit

### **Built-In Tension system**

Option: Each product is sold separately: foot, sleeve, calf, knee and thigh

Warranty: 2 months / factory default

The reduction kits for the legs and arms are for the treatment of acute phase decongestion.

The versatility of the product allows you to meet the specific requirements of your patient. A unique approach with versatile components for a custom solution and repeatable and consistent compression.

For the measurement chart, ask for your form or go to medicanada.ca to download it.



The CircAid Reduction Kit lower leg and armsleeve are compatible and a great addition to the new reduction Kit Glove and Toe Cap.

Learn more on page 58!





Each system includes: 1 Built-in Pressure guide card (BPS) 1 Black Comfort™ Coverup

Lower Leg Upper Leg Full Leg


1 Black CircAid Compressive Undersock 1 Comfort™ PAC Band 1 Juxtafit® Essential Lower leg 1 Comfort™ thigh leg liner(without foot)pair 1 Juxtafit® Essential knee and thigh The full leg is a combination of the lower leg garment and the upper leg with knee garment.

### Juxtafit® Essentials (Lymphedema)

### **Built-In** Pressure **system**

### Warranty: 6 months / manufacturing defect

Juxtafit® Essentials is a line of easy, cost effective ready-to-wear compression leggings for the treatment of lymphedema. Easy to put on and take off, the Juxtafit® Essentials line offers an expanded size range to include 8 sizes with 3 full calf options. Utilizing a comfortable and breathable design, these leggings also feature the same Juxtafit® Premium benefits of anti-bacterial and anti-odor treatments.



Length

Calf Circumference		S	М	<b>M</b> large calf	L	<b>L</b> large calf	XL	XL large calf	XXL
5 cm above ankle	В	15 - 23	20 - 28	20 - 28	25 - 33	25 - 33	30 - 38	30 - 38	33 - 41
15 cm above ankle	<b>B</b>	20 - 28	25 - 33	33 - 40	30 - 38	38 - 46	33 - 41	41 - 51	43 - 53
25 cm above ankle	G	30 - 38	36 - 43	43 - 50	41 - 48	48 - 56	41 - 51	51 - 61	53 - 64
Length	BD		ankle	to the bend o	of the knee le	ngth Short: - 3	6 cm Long: +	36 cm	
Thigh / Knee Circumference		XS	S	М	L	XL			
5 cm below knee	D	27 - 35	31 - 41	35 - 46	40 - 52	44 - 58			
5 cm above knee	<b>(E)</b>	37 - 48	42 - 55	50 - 65	58 - 76	66 - 87			

Knee to groin length Short: - 35 cm Long: + 35 cm





STANDARD

сиѕтом

Standard leg Custom Made leg

The package includes:

 ${\tt 1\,Black\,CircAid\,Compressive\,Undersock}$ 

1 Comfort™ PAC Band

- 1 Juxtafit® Premium Knee, Thigh and/or Lower leg
- 1 Built-in Pressure guide card (BPS)
- 1 Black Comfort™ Coverup

# CIRCAID JUXTA FIT LEGGING BP5" Guide Card 30 40 CUSTOM 40-48cm Mate Course V USE SAME LEVEL FOR ALL BANCS

### Juxtafit® Premium (Lymphedema)


Classes: 20-30, 30-40 and 40-50 mmHg

Custom Made: Choice of Ankle, Lower leg, knee, Thigh or Complete leg

Warranty: 12 months / manufacturing defect

Juxtafit® Premium is the leading inelastic, instantly adjustable compression garment for the management of lymphedema. The easy design enables patients to live independent lives while alleviating the struggle of using bandaging or compression stockings. Made from Breath-O-Prene, the Juxtafit® is soft and comfortable and can be worn day and night. The rigid construction of the Juxtafit® Premium enables it to bridge the gaps associated with severe edema helping to better manage the lymphedema. In addition Juxtafit® Premium has a life span of up to 2 years if cared for appropriately. This is a significant value over other compression garments that may need to be replaced every 6 months.

Calf Standard Circumference		S	М	<b>M</b> large calf	L	<b>L</b> large calf	XL	<b>XL</b> large calf	XXL
Ankle	B	15 - 23	20 - 28	20 - 28	25 - 33	25 - 33	30 - 38	30 - 38	33 - 41
15 cm above ankle	<b>B</b>	20 - 28	25 - 33	33 - 40	30 - 38	38 - 46	33 - 41	41 - 51	43 - 53
25 cm above ankle	G	30 - 38	36 - 43	43 - 50	41 - 48	48 - 56	41 - 51	51 - 61	53 - 64
Foot arch	A1	19 - 24	24 - 29	N/A	29 - 34	N/A	N/A	N/A	N/A
Length	B-D		ankle	to the bend o	of the knee le	ngth short: - 3	6 cm Long: +	36 cm	





The package includes:

- 1 Juxtafit® sleeve
- 1 Built-in Pressure guide card (BPS)
- 1 Comfort™ Coverup (choice of Beige, Black and Pink)



CUSTOM MADE

### **Juxtafit® Lymphedema Essentials Armsleeve**

Classes: 10-20-30-40 mmHg

Option: Comfort arm liner

Warranty: 6 months / manufacturing defect

Tired of struggling with elastic compression sleeves everyday? Then the Juxtafit® armsleeve is the perfect alternative enabling the patient to don the sleeve with a single hand. The Juxtafit® armsleeves provides the compression management needed to maintain a healthy and active lifestyle. The unique profile design ensures perfect Compression in Non-elastic material for a high working pressure. Easy to adjust when the garment is on, the BPS guide card controls the level of compression required.

**Optional:** 



Removable dorsum strap for glove

### **Juxtafit®** Essentials Glove

Warranty: 6 months / manufacturing defect

**Option:** Custom Dorsum Strap



The custom Juxtafit® Glove provides coverage of the hand, wrist and fingers. The garment can adapt to various forms of hands and facilitates the application for the disabled hand.

Option: The Removable Juxtafit® Dorsum strap is a perfect complement to custom Juxtafit® Glove . It can be adjusted to provide more compression on the back of the hand.



## CircAid® by medi

Juxtafit® Gauntlet

Accessories

Juxtafit\*Gauntlet

Available in 3 standard sizes

Comfort<sup>™</sup> Anklet

Accessories

Sold as a pair

Sizes: Standard and Large



Comfort™ PAC Band

Accessories

Designed to be combined with a Comfort™ Anklet or compression stocking

Sizes:: Standard and Large



Comfort™ EZ

Accessories

Width: 3" (7,5 cm)

This product does not contain latex.



Juxtalite® AFW

Juxtafit® AFW

Silvertec treatment

Sizes: Small and Medium

For severe edema Sizes: Small, Medium and Large



Comfort™ Capri

Accessories

Worn under lower leg CircAid\*

Sizes: A to |



Comfort™ Coverup

Accessories

Styles: Sleeve, Lower leg and Full Leg

Sizes: Small and Large



Reduction Kit Toe Cap & Glove

Accessories

Light to moderate compression

Sizes: S to XX-Large



Comfort™ Liner for arm or leg

Accessories

Protects skin

Styles: Sleeve, Sock with or without foot



**CircAid** Compressive Undersock

Protects skin and offers compression to the foot

Accessories

Style: Sock with foot





### **Juxtafit® Essentials Gauntlet**

### Style: Left or Right

The Juxtafit® Gauntlet provides hand coverage from the wrist to the base of the fingers. The gauntlet can be used in combination with a compressive glove. The garment is adjustable and available in five standard sizes, in right or left.

### **Comfort**<sup>™</sup> **Anklet**

Sold as a pair

Sizes: Standard and Large

This product does not contain latex.

The Compression ankle should range around (8-12 mmHg) gives a slight compression on the foot and ankle and can be worn with a lower legging and or CircAid® Comfort™ PAC Band



		Standard	Large
Ankle	B	19 - 27 cm	25 - 50 cm
Foot arch	A1	19 - 36 cm	22 - 41 cm







### **Juxtafit® Premium Interlocking**

#### Sold individually

Sizes: Standard, medium and Large

Adaptive high-quality compression garment for more compression of the foot

The circaid juxtafit premium foot with its individually adjustable foot length and closed heel offers optimum care for the ankle joint and foot. The ankle band provides specific compression of the heel and ankle joint. This foot option is available in three sizes.

### **Comfort**<sup>™</sup> PAC Band

Sizes: Standard and Large

The comfort Power added compression band (PAC<sup>TM</sup>) offers instanly adjustable compression targeting the ankle and arch area of the foot. The Comfort PAC Band<sup>TM</sup> can be trimed to fit almost every foot shape. The lightweight and breathable material fits easily in a pair of shoes, and the comfort PAC Band<sup>TM</sup> is designed to be paired with the comfort compression Anklet or a compression stocking for added compression. Fits Left and Right.

		Standard	Large
Ankle	В	19 - 27 cm	25 - 50 cm
Foot arch	A1	19 - 36 cm	22 - 41 cm



		Standard	Large
Ankle	B	19 - 27 cm	25 - 50 cm
Foot arch	A1	19 - 36 cm	22 - 41 cm







### **Comfort**™ **EZ**

Width: 3"(7.5cm)

Length: 80" (203cm)


#### This product does not contain Latex


The comfort Ez Single-Band-Foot Wrap $^{TM}$ , can be cut to fit any foot size. This wrap gives compression to the foot and ankle. For the difficult foot to fit.

### Juxtalite® Ankle AFW

Sizes: Small and Medium

The ankle-Foot wrap gives ankle and foot coverage. The wrap comes in two sizes and can fit almost any foot. The Breathe-O-Prene material is cool, light-weight, and allows the skin to breathe. The Silvertec<sup>TM</sup> fabric treatment can help eliminate odor and the growth of bacteria. Also, the Juxtalite® ankle-foot wrap is machine washable.









### Juxtafit® Ankle AFW

Sizes: Small, Medium and Large

### For severe edema

The Juxtafit® ankle foot wrap provides ankle and foot coverage. The wrap comes in three standard sizes and can also be custom-made to fit almost any foot shape.

### **Comfort**<sup>™</sup> Capri

### Sizes: A to J


The Comfort™ Capri can be worn with any compression legging for added comfort. Provides support from the waist to just below the knee.

		Small	Medium	Large
Foot arch	A1	19-24 cm	24 - 29 cm	29 - 34 cm
Ankle	B	< 36 cm	< 36 cm	< 41 cm
Length	Z	> 16 cm	> 16 cm	> 16 cm
Lenght	BI	> 9 cm	> 9 cm	> 9 cm



CUSTOM MADE

Circumference							
Waist		Thigh		Waist	Thigh		
A	63 - 79 cm	46 - 62 cm	F	85 - 108 cm	63 - 84 cm		
В	67 - 85 cm	49 - 65 cm	G	90 - 113 cm	67 - 90 cm		
C	72 - 91 cm	52 - 70 cm	Н	94 - 119 cm	70 - 93 cm		
D	76 - 96 cm	55 - 73 cm	1	99 - 125 cm	72 - 96 cm		
E	81 - 102 cm	60-81 cm	J	103 - 130 cm	77 - 102 cm		







### reduction kit™toe cap

Sizes: Small, Medium, Large, X-Large and XX-Large

Light to moderate compression

Patients have the ability to self-manage their lymphedema with ease and confidence with the reduction kit toe cap. Worn with a circaid foot piece, the toe cap provides supplemental decongestive efforts for edema patients.

circumferences in cm.							
small	XX-large						
20-22	22-24	24-26	26-28	28-30			





Trimmable lengths for a perfect custom fit!

### reduction kit ™ glove

Sizes: Small, Medium, Large and X-Large and XX Large

Light to moderate compression

Worn with our inelastic arm & hand wraps, the circaid reduction kit glove effectively improves patient self-management during decongestion. Patients now have the ability to self-manage their lymphedema with confidence and independence.

	ci	rcumferen	ces in cm	
small	medium	large	x-large	xx-large
18-21	20-23	22-25	24-27	26-29









**Up:** Lower leg CoverUp and Full Leg CoverUp

Under: sleeve CoverUp



### Comfort™ CoverUp for leg and sleeve

Sizes: Small and Large

Colours: Beige, Black and Pink

Styles: Armsleeve, Lower Leg and Full Leg

The Comfort™ CoverUp can be worm over the JuxtaLite® or JuxtaFit® garments. Its sleek design prevents the garments from snagging clothing or irritating other parts of the body. With multiple colours available it gives the compression garment a more re ned look.

Note: The comfort  $^{\text{TM}}$  CoverUp does not provide any therapeutic compression.

Circumference	Small	Large
Sleeve CoverUp	26 - 39 cm	40 - 60 cm
Lower Leg CoverUp	34 - 50 cm	51 - 80 cm
Full Leg CoverUp	50 - 71 cm	72 - 108 cm


### **CircAid®** Compressive Undersock

Sizes: Available in standard and large

Styles: Available in 15 - 25 mmHg (moderate compression) and 25 - 35 mmHg (firm compression)

- Applies graduated compression to the foot and ankle
- Graduated compression continues into the leg portion of the undersock
- Soft material for pleasant wearing comfort

	Standard	Large	Size	Strength
сC	26 – 54	41-64	Standard	Moderate
сВ	19-27	25-42	Large	Moderate
cA1	19-32	24-41	Standard	Firm
			Large	Firm




Ordering products has never been so easy!

# medi e-shop

Order online, easily.





- Order standard and custom-made products from the extensive medi range (also by EAN scan).
- Product images and descriptions assist you in the ordering process.
- Integrated size checks and size information make queries and complaints less likely.
- See expected shipment date for your order.
- Carry out follow-up orders for your customers with just a few clicks.
- Discover our digital browsing catalogues.

Visit **eshop.medicanada.ca** today!

medi. I feel better.

Notes	





medi Canada inc 597, rue Duvernay Verchères QC Canada JOL 2R0

service@medicanada.ca medicanada.ca 1 800 361-3153